



TOSOH


-sizing guide

OCTAVE BIO & OCTAVE PRO



Octave BIO & PRO

Catalog No.	Systems	Min. Flow Rate	Max. Flow Rate	Column Volume	Max Capture Throughput
0041100	Octave BIO 12	0.1 mL/min	12 mL/min	<25 mL	170 g/day
0041101	Octave BIO 36	0.1 mL/min	36 mL/min	<72 mL	520 g/day
0041102	Octave BIO 100	1 mL/min	100 mL/min	1-200 mL	1.4 kg/day
0041103	Octave BIO 300	3 mL/min	300 mL/min	5-600 mL	4.3 kg/day
0041200 US 0041201 EU	Octave PRO	100 mL/min	2.5 L/min	0.2-5 L	36 kg/day



TOSOH BIOSCIENCE
**SEPARATION
& PURIFICATION**
CONNECTING MINDS.
TOUCHING LIVES.

For more information:

Contact our #ChromatographyExperts at
www.tosohbioscience.com

Accessories

Catalog No.	Accessories	Octave BIO 12	Octave BIO 36	Octave BIO 100	Octave BIO 300
0041104	1 mm Ultem Valve Block	1x included	1x included		
0041105	1.6 mm Ultem Valve Block			1x included	
0041106	2.4 mm Ultem Valve Block				1x included
0041107	Octave 12 mL/min pump	6x included			
0041108	Octave 36 mL/min pump		6x included		
0041109	Octave 100 mL/min pump			6x included	
0041110	Octave 300 mL/min pump				6x included
0041123	Octave, Tubing, FEP, 1/16" OD, 0.02" ID x 50'	1x included			
0041124	Octave, Tubing, FEP, 1/16" OD, 0.03" ID x 50'		1x included		
0041125	Octave, Tubing, FEP, 1/8" OD, 0.062" ID x 50'			1x included	
0041126	Octave, Tubing, FEP, 1/8" OD, 0.093" ID x 50'				1x included
0041145	Octave BIO 1 mm valve block, PCTFE				
0041146	Octave BIO 1.6 mm valve block, PCTFE				
0041147	Octave BIO 2.4 mm valve block, PCTFE				
0041127	Octave, Tubing, PFA, 1/16" OD, 0.02" ID x 50'				
0041128	Octave, Tubing, PFA, 1/16" OD, 0.03" ID x 50'				
0041129	Octave, Tubing, PFA, 1/8" OD, 0.062" ID x 50'				
0041130	Octave, Tubing, PFA, 1/8" OD, 0.093" ID x 50'				

Additional accessories can be purchased separately to convert an Octave BIO system to a different size.

TOSOH BIOSCIENCE

**SEPARATION
& PURIFICATION**

CONNECTING MINDS.
TOUCHING LIVES.

For more information:

Contact our #ChromatographyExperts at
www.tosohbioscience.com